Instructions for use

Torque Wrench

mis[®]

Description of the torque wrench

The torque wrench with adjustable force is a dental device used to tighten or loosen screws, prosthesis components and implants. It is a precision instrument that can be disassembled and that is supplied unsterile. To ensure that it functions perfectly every time, the torque wrench must be disassembled, disinfected, cleaned, greased and sterilised after each use, according to the instructions for use. Persons who use this device are obliged to read the instructions for use prior to handling. The handling and the use of the product are carried out without direct control from our side and remain under responsibility of the user. The user is liable for any possible damage that could occur. Before each use, in order to guarantee high torque precision, the device must be checked upon its functioning. This instrument is not a measuring device.

Use

By turning the torque adjustment screw, the torque wrench can be set to the desired torque value. To set the torque value correctly, the torque adjustment screw must be turned clockwise to reach the required torque value and set to the exact line marking. Ensure that the line on the handle is in straight alignment with the line on the torque adjustment screw. In order to change from a higher to a lower torque value, one must screw two turns under the desired torque value, then screw clockwise again to the exact line marking. Ratchet mode could be set by turning the torque adjustment screw until the R marking for the torque wrench including Ratchet option. The word 'IN' on the cover (3), shows the position of the wrench that is used for tightening ; by turning the device over, the word 'OUT' is used for loosening screws.

Grease

"Instrument Lubricant" approved USDA H1


Precision of new device

± 3,5 Ncm with total confidence of 95%

Recommendations

This instrument must not be used for any applications other than those listed in the section "Description of the torque wrench" or with equipment that could damage the intended use of the device. The persons in charge for the use and maintenance of this dental instrument should monitor any deterioration of the tightening, ratchet and torque mechanism of the device and, in the event of a defect, return the wrench to the supplier. During assembly, it is essential not to mix the various components belonging to different torque wrenches because the components are not interchangeable. If a component is lost, please return the whole instrument immediately to your retailer for repair. Components cannot be sold separately. Do not store the wrench with the spring compressed but with the torque set to its minimum. This device must not be sterilised in the packaging provided by the manufacturer


Individual parts of the torque wrench


Cleaning the torque wrench

When used in situations that do lead to operative residues (blood, secretions, tissue remnants), the torque wrench must be disassembled completely (fig. 2 to 4) and placed in a suitable bath of disinfection in accordance with the recommendations of the manufacturer. This operation facilitates cleaning because dry residues cause corrosion. After cleaning, thoroughly rinse the parts with water and use a nylon brush to rub internal and external surfaces of the various parts of the torque wrench. During the cleaning process, avoid all contact between each part of the torque wrench.

Disassembling the torque wrench

Unscrew completely the torque adjustment screw (7) and extract the spring/ stop (6). If necessary, use the hexagon shank bit (8) to extract the spring. Be careful not to damage the hexagon shank bit. Do not separate the stop from the spring (6).


With the hexagon shank bit (8) of the torque adjustment screw (7), unscrew the screw (5), while maintaining a light pressure, to remove the cover (3).


The ratchet head is composed of two removable elements: the ratchet (2) and the ratchet wheel (1). The screw (5) is fixed inside the ratchet head.

Assembling the torque wrench

To assemble the torque wrench, insert the next two elements in the following sequence: first the ratchet wheel (1) and then the ratchet (2).

Moderately lubricate the functional contact area between the teeth of the ratchet wheel (1) and the pivot point of the ratchet (2) as shown. Remove any trace of the lubricant from the outer surface of the torque wrench because excessive lubricant leaves traces on the instrument's surface during sterilisation. Use only the instrument lubricant delivered with the torque wrench.

Once the parts (1) and (2) are in place, put the cover (3) on the ratchet head (4) and adjust. Turn the torque wrench while holding the cover. Fix the screw (5) tightly with the hexagon shank bit (7) of the torque adjustment screw.

The assembly of the torque wrench is complete as soon as the spring/stop is in place and the torque adjustment screwed back in the handle. Check the function by activating the ratchet wheel. If there is a regular ratchet noise, the instrument operates correctly. Before use, always check the assembly and the correct function of the torque wrench.

Use the adjustment key (provided) to obtain the higher torque values more easily.

Sterilization

The instrument must undergo a sterilization with steam at 134 °C/ 273°F during 18 minutes. Before sterilization, the torque wrench must be completely assembled (5-8). Sterilise the key according to cycles of sterilization recommended by the manufacturer of the autoclave. We recommend the use of devices equipped with a vacuum pump (type B) to decrease the risk of formation of air pockets. This recommendation is particularly important for hollow instruments and guarantees a perfect drying. We advise against the use of a hot air steriliser because it can lead to ageing of the spring and subsequently bring about a change of the torque value.

Key to codes used

Attention, see instructions for use

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All rights reserved. No part of this publication may be reproduced, transcribed, stored in an electronic retrieval system, translated into any language or computer language, or be transmitted in any form whatsoever, without the prior written consent of the publisher. Warning: Only a licensed dentist should use these products.
MIS Implants Technologies Ltd.
MIS Implants Technologies GmbH