


MIS System Guide | 2010

**MIS<sup>®</sup>**  
**User Manual**



© MIS Corporation. All Rights Reserved.


Published by MIS, which reserves the right to modify the products described in this manual as well as to revise this publication at any time and without informing any person of such revision or change. All rights reserved. No part of this publication may be reproduced, transcribed, stored in an electronic retrieval system, translated to any language or computer language, or transmitted in any form whatsoever without the written consent of the publisher. Questions, comments or requests will be addressed promptly by contacting MIS specialists directly through our e-mailing address: [service@mis-implants.com](mailto:service@mis-implants.com). MIS' website can be accessed at [www.mis-implants.com](http://www.mis-implants.com). This online site highlights current products and reflects all new discoveries and developments.

Note: This User Manual is for educational use only.

MIS Implants Technologies Ltd. produces an oral implant system that includes self-tapping dental implants in a wide range of sizes. These implants, currently being used successfully worldwide, provide a solution to partial and full edentulism. All implants go through gamma-irradiation and are supplied sterile in specially designed tubes. The implants are self-tapping with threads designed to provide secure primary fixation and favorable distribution of the loading forces. The MIS Implants system is manufactured from titanium alloys. MIS Implants Technologies Ltd. is vigilant about maintaining the high quality of its products as well as developing new products in the fields of dental implants, abutments and surgical kit components.

---

MIS's Quality System complies with international quality standards: **ISO 13485:2003** - Quality Management System for Medical Devices, **ISO 9001: 2008** – Quality Management System and CE Directive for Medical Devices **93/42/EEC**. MIS's products are cleared for marketing in the USA and CE approved.





**Overview.**  
MIS concept

6. Introduction
7. Raw Material
8. Finite Elements Test
10. Implant Treatment
11. Implant Surface
13. Roughness Measurements
14. XPS Analysis
15. In Vitro Test
16. Histology

## Overview

# Introduction

MIS is a dynamic, high-tech research and production company. It develops and manufactures a comprehensive range of dental implants that provide long-lasting successful solutions to partial and full edentulism. The MIS Implants system combines several advantageous elements in order to achieve good primary stability and successful osseointegration. The qualities that contribute to successful osseointegration have been determined by various researchers to include: choice of raw material, macrostructure, microstructure and surface cleaning. In this chapter we present these factors and others that have been considered and/or included in the implants manufacturing process. MIS upholds its high standards, relevant to these factors, through comprehensive examinations and tests.

The implant surface is tested and has satisfied or exceeded all requirements for implants. The tests are carried out in some of the world's best-known research institutions. The main tests that MIS implants undergo are:

- Mechanical test
- XPS Analysis
- Roughness
- Surface analysis
- SEM
- Cytotoxicity
- Sterility
- Torque removal value
- Histology

7.

## Overview Raw Material

All MIS implants are made of biocompatible medical grade Titanium alloy.

Titanium and its alloys are highly successful materials for the fabrication of dental implants because of their favorable combination of properties such as; low specific weight, high strength to weight ratio, high modulus of elasticity, very high corrosion resistance and excellent general biocompatibility. The excellent biocompatibility and osseointegration capabilities of titanium are related to a variety of favorable properties of the material and its surface, including the following:

- A dense, highly resistant passive oxide film that protects the underlying metal from (further) oxidation and corrosion.

- A very low dissolution rate of the oxide film and an extremely low concentration of charged titanium corrosion products.

Biocompatibility is a complex property that involves physical, chemical, biological, medical and design aspects. It is important to note that the corrosion resistance of the titanium alloy is directly related to the surface condition before corrosion evaluation, and is a determining factor in the electrochemical behavior of this alloy. A study of titanium and its alloys, demonstrated that the biocompatibility of metal implants was related to their electrochemical characterization. The results showed that the presence of cells on metals led to an increase in

the impedance and polarization resistance ( $R_p$ ) values of the metals. SEM micrographs, and an equivalent circuit confirmed this behavior.


Modification of the surface by thermal treatment improves corrosion resistance, as a protective titanium layer is formed on the surface. Another quality that makes titanium uniquely suited to implants is that the titanium oxide layer is stable over a very large range of physiologic pH, and therefore there is no release of free ions into the tissue. Titanium alloys increase the mechanical strength by 30%-40% over pure titanium.



◀ Titanium alloys increase the mechanical strength by 30%-40% over pure titanium.

Overview  
**Finite  
Elements Test**


The finite elements method (FEM) provides additional information about implants. Implants analysis by FEM is a test that is commonly used for the localization and prediction of a mechanical system failures. A structure analysis enables the determination of effects such as deformations, strains and stresses which are caused by applied structure loads. The safety factor is calculated using the maximum stress failure theory for ductile materials. The stress limit is specified by the tensile yield strength of the material.




The benefits of FEM include increased accuracy, enhanced design and better insight into critical design parameters. The results of these tests show that the greatest stress occurs in the neck rather than the apex of the implant.



MF4



MF7



## Overview **Implant Treatment**



Structure (Raw Material)



Roughness  
(Macro & Micro)

Titanium's rough surface provides many points for mechanical anchorage and increases the envelope surface by two to three times.



Subtraction

Advantage: no degradation and separation of the coating.



Sandblast

The combination of the two methods induces macro and microstructure that is optimal for cell adhesion.



Acid Etching




MIS Surface Treatment


11.

## Overview Implant Surface

The contact surface of a dental implant is a very important aspect of the implant's long-term success. This phenomenon is known as osseointegration. The opposing reaction of the bone to the implant is determined by morphology, chemical composition, load and the characteristics of the implant surface. It is achieved by topography and morphology, nano surface elements and bone reaction tests.



SEM image of the bone-to-implant interface



SEM image showing an osteoblast

At MIS all of the following aspects have been considered:

### **Macrostructure**

The geometric design of the body and thread profile of the implant act to increase the primary stability and to distribute the axial forces of occlusion.

### **Microstructure**


The surface roughness and microgeometry of titanium are achieved by sandblasting and acid-etching. This process increases the surface envelope of the implant. Blasted surfaces demonstrate more bone to implant surface contact when compared with machined surfaces.

\*This advanced geometric design and surface morphology are the factors that create successful implants.


### **Surface cleaning**

The implant surface consists of titanium oxide (Ti O<sub>2</sub>), which can reach a thickness

of 100 microns, and is in turn covered with an additional one-molecule-thick layer of contaminants containing sodium, silicon, magnesium, fluorine, etc. The cleaning process is multi- staged in order to ensure the quality of the product. The success of the implant is determined its geometric shape, by its microscopic texture and its surface quality.



x2000




x5000

13.

## Overview

# Roughness Measurements


Roughness or rugosity is a measure of the small-scale variations in the height of a physical surface. A rough texture allows the trapping of different elements and increase the envelope surface. A high tendency for high roughness improves mechanical anchorage and increases the envelope surface by two to three times.



**AFM Test**  
MIS Surface topography

### Roughness Measurement

The arithmetic average of the deviation  $R_a$  is the most commonly used measure for surface roughness. The microgeometry of MIS implants meets the roughness recommended in the international literature.



**Instrument:**  
Parthometer M1  
(MAHR)

$R_a$  2.2 $\mu$ m  
 $R_z$  14.6 $\mu$ m  
 $R_{max}$  15 $\mu$ m  
 $L_t$  5.6mm  
 $L_c$  0.8mm  
 $P_c$  (0.5-0.5) 165/c

## Overview XPS Analysis

For surface analysis the samples were irradiated with monochromatic X-rays. Survey spectra were recorded with a pass energy of 100 eV, allowing the surface chemical composition to be determined. The atomic concentrations were calculated using elemental sensitivity factors,

without applying any standardization procedure. The core level binding energies of the different peaks were normalized by setting the binding energy for the Cls to 284.5 eV. For each of the screws, one flat edge area, marked as #1, was analyzed in the as-received state only.

**Table 1** XPS Atomic Concentrations (%) for LOT 83

Implants	C	Ti	O	N	Ca	Si	S	Cl	Na	Al	Cu	Mg	P	V
MF7-13375, thread area 1	24.31	19.45	52.71	-	-	-	0.12	-	3.08	-	-	-	-	0.33
MF7-13375, thread area 2	27.78	16.11	51.31	0.21	-	0.05	-	0.78	3.60	-	-	-	-	0.16

**Table 2** Ti Oxide Thickness (nm) for LOT 83

Samples	Oxide Thickness (nm)
MF7-13375, thread area 1	6.5

**Table 1** XPS Atomic Concentrations (%) for LOT 41903


Implants	C	Ti	O	N	Ca	Si	S	Cl	Na	Al	Cu	Mg	F	P	ZN
MF7-13375, thread area 1	25.71	17.60	52.22	1.99	-	-	-	-	1.96	-	0.51	-	-	-	-
MF7-13375, thread area 2	28.66	14.88	52.62	0.84	-	-	-	--	2.12	-	0.87	-	-	-	-

\*Instrument: VG Scientific Sigma Probe  
X-Ray Source: Monochromatic Al K $\alpha$ , 1486.6eV  
X-Ray Beam Size: 400  $\mu$ m

### Interpretation of Results

No strange elements or traces thereof were identified. This means that the surface treatment (etching) did not leave any undesired effects. The implants were classified as normal based on the C/Ti ratio.

## Absorption of Serum Protein to Modified Titanium Surfaces



**M.N. Sela, L.Badihi, G.Rosen  
D.Kohavi and D. Steinberg**

The use of Titanium (Ti) implants is an innovative clinical procedure in dentistry. The absorption of biological molecules to the implant's surface triggers a sequence of events that can determine the outcome of this procedure. Clinical data suggests that modified Ti surfaces play an important role in the success or failure of the implant.

**Objective:** the purpose of this study was to investigate the interaction between Ti implants with different surface properties and serum proteins, in order to find the optimal implant surfaces that may improve the osseointegration process and implant intake.

**Materials & Methods:** Ti disks (6mm in diameter) with two types of surface modifications were compared: Machined and Sandblast together with Acid-Etched. The disks were coated with mixtures of Human Serum Albumin conjugated with fluorescein (HAS-FITC). Following incubation, the coat was removed from the disks by SDS.

A Confocal Scanning Laser Microscope was used to visualize and measure the HAS-FITC coat and the degree of protein removal from the Ti surfaces.

**Results:** The Confocal Microscope images revealed a significantly higher amount of HAS-FITC coat on the rough disks, as compared with the machined disks. Furthermore, under similar experimental conditions, less HAS-FITC could be removed from the rough disks than from the machined disks.




**Conclusions:** Absorption of albumin to the rough treated Ti surface is both qualitatively and quantitatively far more intense, as compared with the machined surfaces. Further studies of the chemical and physical characterization of the modified Ti surfaces are underway. Moreover, additional serum proteins, as well as oral microorganisms, are being examined for their interactions with the modified Ti surfaces.

16.

Overview  
**Histology**

MIS implants undergo routine testing, including histology and microsimy tests, in order to evaluate the amount and quality of bone integrated on the implant as well as other indications for osseointegration. The osseointegration's success is determined through the many contact areas of the implant surface. The following report characterizes the MIS implant surface as indicated by the various test results.

**Light Microscopy**





## Materials and Methods

### Micro CT

The sample was analyzed by micro CT before histological treatment was analyzed by micro CT. X-ray micro-computed tomography (SkyScan 1072, Belgium). The X-ray source was operated at 100 kV and a current of 98  $\mu$ A. The sample was rotated 180° with a rotation step of 0.90°, an acquisition time of 5.6 s per scan and a pixel size of 11.8  $\mu$ m.

Three-dimensional reconstructions were then performed with the 3D Creator SkyScan software.

### Histology

#### Embedding

Samples were fixed in a 10% formalin solution, dehydrated in ascending graded ethanol (70, 80, 90, 95, 100%) and pure acetone, and then impregnated and embedded in PMMA (methyl methacrylate).

The Polymerization was obtained with a temperature increase from 20°C to 80°C.

## Results and Discussion

The histology and microscopy results show the following osseointegration contact areas between the bone and implant surfaces:

### Light Microscopy

Full osteointegration was observed. Bone contact was achieved both in the front of the wire and at the top of the implant. Bone ingrowth with bone trabeculae was also observed.


### Micro CT

Bone osseointegration was observed along all the 3D axes. Bone trabeculae are directly in contact with the wires.

### Conclusion

In conclusion, it appears that larger bone ingrowth and bone contact was observed in the implants.

## X-ray microtophography







**Implants.**

- 20. Implant System
- 22. UNO
- 30. BIOCOM
- 38. SEVEN

20.

# Implants

Implant System

The MIS leading implants:  
UNO, BIOCOM and SEVEN



MIS has developed a range of unique implants and tools to assist in the simplification of the implantation process and to ensure both efficiency and success while minimizing risks.

The innovative design of our implants, combined with our simple and fast insertion procedures, provides an easy to use system that ensures successful results.

The wide range of MIS implant lines provides a variety of clinical solutions for the reconstruction of a single tooth, screw-retained or cemented fixed bridges, and overdentures. Furthermore, MIS

implants can be used in any surgical and bone augmentation procedures, from the simplest to the most intricate. MIS implants are made of high-quality materials under very strict quality control procedures with a 100% inspection rate.

All MIS implants are made of biocompatible medical-grade titanium, and their surfaces undergo dual-acid-etching procedures.



**Note:** The selection of MIS implants and procedures is a recommendation only, and cannot replace the judgment and the experience of the surgeon. This recommendation should be used as a guide only.

## UNO Introduction

Each of these unique implants is specifically engineered for use in narrow ridges and tight spaces. The insertion of the UNO implant is a quick and simple one-stage procedure. Due to their innovative geometries and advanced surface morphology, these unique implants offer high primary stability. These versatile implants can be used to restore single crowns and anterior cemented bridges.


UNO NARROW ▶



◀ UNO ONE PIECE




UNO One Piece  
Fixture - Technical Info



- 1 Abutment with squared crown
- 2 Deep threads 0.6mm
- 3 Dual thread of 2mm
- 4 Cutting Apex
- 5 Two spiral channels
- 6 Implant length
- 7 Surface: (A) acid etched (B) sandblasted and acid etched
- 8 Maximum angled position 10°


# UNO NARROW

Fixture - Technical Info



- 1 External Hexagon 2.7 mm
- 2 Deep threads 0.6mm
- 3 Cutting Apex
- 4 Two spiral channels
- 5 Implant length
- 6 Dual thread of 2mm
- 7 Surface- sandblasted and acid etched





## UNO Features

### Features

- A mono block screw type implant. Especially designed for narrow ridges and tight spaces.
- The geometric is encompassed in the unique design.
- Bone-to-implant contact maximized.
- Enables a fast, thread-form free insertion and initial stability.

### Advantages

#### Simple

MIS UNO's specially designed tools and simple procedure ensure a worry-free restoration for the clinician.

#### Easy

An innovative design with increased insertion speed making the MIS UNO an easy implant to insert.

#### Stability

The MIS UNO design ensures maximum strength and stability for the implant and restorative parts.

#### Versatility

The MIS UNO is indicated for use in narrow ridges and tight places, such as maxillary lateral and mandibular incisors. The MIS UNO's versatility enables the clinician to use the implant for single-tooth, partial-denture and over-denture restorations.

#### Long-lasting









Due to innovative geometry and advanced surface treatment, the MIS UNO provides high initial stability and a long-lasting restorative result.

#### Surface

The surface roughness and microgeometry of Titanium Alloy Ti 6Al-4V ELI are achieved by sand blasting and acid etching. Blasted surfaces demonstrate more bone to implant surface contact compared with machined surfaces.

## UNO Implant Range

One piece -----

Length	10mm	11.50mm	13mm	16mm
Type	MO1-10300	MO1-11300	MO1-13300	MO1-16300
<b>Ø3mm</b> One Piece Implant				
	MO1-10350	MO1-11350	MO1-13350	MO1-16350
<b>Ø3.50mm</b> One Piece Implant				

Narrow -----

	MO2-10300	MO2-11300	MO2-13300	MO2-16300
<b>Ø3mm</b> Narrow Two Piece Implant				
UNO Ball Anchor Screws	MB-NR005  0.50mm	MB-NR015  1.50mm	MB-NR025  2.50mm	


**UNO**  
 Ø3mm / Ø3.5mm  
 One Piece Procedure

**Ø3mm**

<b>Drill Speed (RPM)</b>	1200- 1500	900- 1200		700- 900	15-20
<b>Diameter</b>	Ø1.90	Ø2	Ø2	Ø2.40	Ø3


**Ø3.50mm**

<b>Drill Speed (RPM)</b>	1200- 1500	900- 1200	700- 900	400- 700	15-20	
<b>Diameter</b>	Ø1.90	Ø2	Ø2	Ø2.40	Ø3	Ø3.50



UNO  
Ø3mm Narrow  
Procedure

Ø3mm	1200- 1500	900- 1200	700- 900	15-25	
Drill Speed (RPM)					
Diameter	Ø1.90	Ø2	Ø2	Ø2.40	Ø3



Any procedure recommended by MIS cannot replace the judgment and the experience of the surgeon.


30.

## BiOCOM Introduction

MIS BIOCOCOM implants are titanium cylinder screw type implants with an internal hexagon connection. They are designed for both two-stage or single-stage procedures. The implant is self-tapping and has a unique wide thread design as well as tapered threads at the apical part.



**BiOCOM**  
Fixture - Technical Info





## BiOCOM Features

### Features

- The BiOCOM implants are self tapping with threads that are designed to provide secure primary fixation.
- A wide range of restoration parts.
- Suitable for two-stage implant procedures.
- Available in 3.30mm, 3.75mm, 4.20mm, 5mm and 6mm diameters and lengths of 6mm, 8mm, 10mm, 11.50mm, 13mm and 16mm (see page 25).

### Surface



























The surface roughness and microgeometry of Titanium Alloy Ti 6Al-4V ELI are achieved by sand blasting and acid etching. Blasted surfaces demonstrate more bone-to-implant surface contact compared with machined surfaces.

### Self-Tapping

The three cutting flutes are designed to engage the bone immediately during placement and ensure multi-directional locking. The tapping head cuts into the bone with far less friction due to the relief design of the cutting edge.



## Implant Range

Length	6mm	8mm	10mm	11.50mm	13mm	16mm
<b>Ø3.30mm</b> Screw type implant narrow platform			MF4-10330 	MF4-11330 	MF4-13330 	MF4-16330 
<b>Ø3.75mm</b> Screw type implant standard platform		MF4-08375 	MF4-10375 	MF4-11375 	MF4-13375 	MF4-16375 
<b>Ø4.20mm</b> Screw type implant standard platform	MF4-06420 	MF4-08420 	MF4-10420 	MF4-11420 	MF4-13420 	MF4-16420 
<b>Ø5mm</b> Screw type implant wide platform	MF4-06500 	MF4-08500 	MF4-10500 	MF4-11500 	MF4-13500 	MF4-16500 
<b>Ø6mm</b> Screw type implant wide platform	MF4-06600 	MF4-08600 	MF4-10600 	MF4-11600 	MF4-13600 	

\* All implants include the cover screw.

## BiOCOM

### Ø3.30mm / Ø3.75mm Procedure

\* Recommended insertion torque is: 35-60 Ncm.

#### Ø3.30mm

**Drill Speed (RPM)**

1200-1500    900-1200    500-700    200-500    15-25

**Diameter**

Ø1.90    Ø2    Ø2    Ø2.80    Ø2.80    Ø3.30    Ø3.30



1

**Countersink**  
(MT-GDN33)  
for bone type 1&2


#### Ø3.75mm

**Drill Speed (RPM)**

1200-1500    900-1200    500-700    400-700    200-500    15-25

**Diameter**

Ø1.90    Ø2    Ø2    Ø2.80    Ø3.20    Ø3.20    Ø3.75    Ø3.75



1

**Countersink**  
(MT-GDN33)  
for bone type 1&2



Any procedure recommended  
by MIS cannot replace the  
judgment and the experience  
of the surgeon.

## BiOCOM

Ø4.20mm / Ø5mm  
Procedure

\* Recommended insertion torque is: 35-60 Ncm.

### Ø4.20mm

<b>Drill Speed (RPM)</b>	1200-1500	900-1200		500-700	400-700	400-600		200-500	15-25
<b>Diameter</b>	Ø1.90	Ø2	Ø2	Ø2.80	Ø3.20	Ø3.80	Ø3.80	Ø4.20	Ø4.20



Countersink  
(MT-GDN33)  
for bone type 1&2

### Ø5mm

<b>Drill Speed (RPM)</b>	1200-1500	900-1200		500-700	400-700	400-600	400-600		200-500	15-25
<b>Diameter</b>	Ø1.90	Ø2	Ø2	Ø2.80	Ø3.20	Ø3.80	Ø4.50	Ø4.50	Ø5	Ø5



Countersink  
(MT-GDN50)  
for bone type 1&2

# BiOCOM

## Ø6mm Procedure


\* Recommended insertion torque is: 35-60 Ncm.

### Ø6mm

Drill Speed (RPM)	1200- 1500	900- 1200	500- 700	400- 700	400- 600	400- 600	300- 500	300- 500	200- 500	15-25		
Diameter	Ø1.90	Ø2	Ø2	Ø2.80	Ø3.20	Ø3.80	Ø4.50	Ø5	Ø5.50	Ø5.50	Ø6	Ø6

**1**  
Countersink  
(MT-GDN50)  
for bone type 1&2





38.


## SEVEN Introduction

The MIS self-tapping SEVEN implants are specially designed for implantation in a wide range of bone types and bone augmentation procedures. Their geometric design includes dual threads, three spiral channels stemming from the apex, micro rings on the implant neck and a changing thread thickness along the implant. All MIS SEVEN implants are supplied with a single use final drill for reducing the heat produced during drilling, which results in an improved osseointegration.




# SEVEN WIDE STANDARD

Fixture - Technical Info



# SEVEN NARROW

Fixture - Technical Info



- 1 **Narrow:** 2.10mm Int. hex. connection
- 2 Conical body
- 3 Surface - sandblasted+acid etched
- 4 Two spiral channels
- 5 Domed apex
- 6 **Narrow:** Dual thread of 2mm
- 7 Micro-rings (0.1x0.3mm)



41.

## SEVEN Features

### Features

- The SEVEN implant is designed to suit a wide range of bone types and bone augmentation procedures.
- A specially designed drill ensures short and safe drilling procedures.
- A double thread of 2.40mm increases the implant's insertion speed.
- It has self-tapping capability.
- It has three spiral channels for improved integration.
- The micro-rings (0.1x0.3mm) on the implant's neck reduce the shear stress in the crest zone.
- The changes in the thickness of the thread (0.15-0.4 mm) improve bone compression.
- The SEVEN implants are available in 3.30mm, 3.75mm, 4.20mm, 5mm and 6mm diameters and 6mm, 8mm, 10mm, 11.50mm, 13mm and 16mm lengths.

### Successful

The SEVEN implant has a high success rate as a result of its advanced geometric design and new surface morphology.

### Forgiving

SEVEN is designed for implantation in a wide range of bone types and bone augmentation procedures.

### Simple

A specially designed final drill is supplied with every implant, allowing a short and safe drilling procedure.

### Easy

Increased insertion speed is provided by a dual thread of 2.40mm, combined with a self-drilling capability.

### Primary Stability

The thread thickness changes from the apex to the neck with the same pitch, improving the compression of the bone during insertion. Micro-rings on the implant's neck provide better initial stability by improving the interfacial shear strength at the crest zone.

### Minimal Bone Resorption



























The surface roughness over the entire body, the unique surface morphology, together with the micro-rings at the implants neck, prevents bone resorption at the implant's neck.

### Self-Tapping

SEVEN cuts its own threads during implantation, minimizing friction-generated heat. Three spiral channels running the length of the fixture fill in with bone chips during implantation to improve integration.

# SEVEN Implant Range

42.

Length	6mm	8mm	10mm	11.50mm	13mm	16mm
<b>Ø3.30mm</b> Screw type implant narrow platform			MF7-10330 	MF7-11330 	MF7-13330 	MF7-16330 
<b>Ø3.75mm</b> Screw type implant standard platform		MF7-08375 	MF7-10375 	MF7-11375 	MF7-13375 	MF7-16375 
<b>Ø4.20mm</b> Screw type implant standard platform	MF7-06420 	MF7-08420 	MF7-10420 	MF7-11420 	MF7-13420 	MF7-16420 
<b>Ø5mm</b> Screw type implant wide platform	MF7-06500 	MF7-08500 	MF7-10500 	MF7-11500 	MF7-13500 	MF7-16500 
<b>Ø6mm</b> Screw type implant wide platform	MF7-06600 	MF7-08600 	MF7-10600 	MF7-11600 	MF7-13600 	

\* All implants include the cover screw and final drill.

43.

**SEVEN**  
 Ø3.30mm / Ø3.75mm  
 Procedure

\* Recommended insertion torque is: 35-60 Ncm.


**Ø3.30mm**

<b>Drill Speed</b> (RPM)	1200-1500	900-1200	500-700	200-400	15-25
<b>Diameter</b>	Ø1.90	Ø2	Ø2.40	Ø2.20 Ø3.20	Ø3.30



**Ø3.75mm**

<b>Drill Speed</b> (RPM)	1200-1500	900-1200	500-700	200-400	300-500	15-25
<b>Diameter</b>	Ø1.90	Ø2	Ø2.80	Ø2.80 Ø3.60	Ø2.80 Ø3.60	Ø3.75



**1**  
**Final drill**  
 For bone type 1&2

**2**  
**Countersink**  
 (MT-GDN3)  
 For bone type 3&4

**⚠** Any procedure recommended by MIS cannot replace the judgment and the experience of the surgeon.

## SEVEN

Ø4.20mm / Ø5mm  
Procedure

\* Recommended insertion torque is: 35-60 Ncm.

### Ø4.20mm

**Drill Speed** (RPM)

1200-1500    900-1200    500-700    400-700    200-400    200-500    15-25

**Diameter**


Ø1.90    Ø2    Ø2    Ø2.80    Ø3.20    Ø3.20    Ø3.30    Ø3.30    Ø4.10    Ø4.20    Ø4.20

1

**Final drill**  
For bone type 1&2

2

**Countersink**  
(MT-GDN33)  
For bone type 3&4



### Ø5mm

**Drill Speed** (RPM)

1200-1500    900-1200    500-700    400-700    400-600    200-400    200-500    15-25

**Diameter**


Ø1.90    Ø2    Ø2    Ø2.80    Ø3.20    Ø3.80    Ø4.10    Ø4.10    Ø5    Ø5    Ø5

1

**Final drill**  
For bone type 1&2

2

**Countersink**  
(MT-GDN50)  
For bone type 3&4



# SEVEN

## Ø6mm Procedure

\* Recommended insertion torque is: 35-60 Ncm.

### Ø6mm

Drill Speed (RPM)	1200- 1500	900- 1200	500- 700	400- 700	400- 600	400- 600	300- 500	200- 400	200- 500	15-25	
Diameter	Ø1.90	Ø2	Ø2	Ø2.80	Ø3.20	Ø3.80	Ø4.50	Ø5	Ø5.10 Ø5.90	Ø6	Ø6

**1**  
Final drill  
For bone type 1&2

**2**  
Countersink  
(MT-GDN50)  
For bone type 3&4





# **Surgical Procedures.**

For MIS Implants

- 48. Indications & Contraindications
- 50. Step by Step Protocol

## Surgical Procedures

# Indications and Contraindications



### Indications

Adequate bone is needed to support the implant with width and height being the primary dimensions of concern. The obligatory bone dimension for implant placement in a chosen site may be extracted for adequate radiological techniques used in implant dentistry.

In addition a very careful evaluation has to be made as to the location of vital blood vessels, nerves, maxillary sinus, soft tissue spaces, and their relation to the site planned for implant placement.



### Contraindications

The contraindications customary in oral surgery with other implant materials should be observed. These include patients taking corticosteroids, anticoagulants or anticonvulsant and those receiving radiation or other immunosuppressive therapy. Lactating or pregnant women are not candidates, nor are patients with abnormal laboratory values for BUN,

creatinine or serum calcium. Patients with diabetes or cardiovascular disease are contraindicated. Hypertension above 110/170 mmHg, osteoporotic crush fractures, respiratory disease, thyroid or parathyroid disease should be excluded from treatment. Patients with diagnosed malignancy in the past five years and those with nodular enlargements, tenderness or unexplained lumps or masses of the head or neck should not be treated. Implanting procedures should not be performed on persons with active osteolytic, inflammatory or infectious processes in the implantation site.

### The following outline lists the contraindications:

- Debilitating or uncontrolled disease.
- Pregnancy, Hemophilia, Granulocytopenia or other bleeding problems, steroid use, Prophylactic antibiotics, Brittle diabetes, Ehler-Danlos syndrome.
- Osteoradionecrosis, renal failure, Organ transplantation anticoagulation therapy.
- Unexplained hypersensitivity, Fibrous dysplasia, Regional enteritis.





## Other Contraindications

Conditions, diseases, or treatments that severely compromise healing, including radiation therapy.

Poor patient motivation.

Psychiatric disorders that interfere with patient understanding and compliance with the necessary procedure.

Unrealistic patient expectations.

Unattainable prosthodontic reconstruction.

Inability of patient to manage oral hygiene.

Patient hypersensitivity to specific components of the implants.

infections in susceptible individuals, including those with body part replacement. Existing natural dentition may be compromised by improper implant placement.

The following outline lists organ systems with corresponding pathophysiological problems that may influence risks:

- a. Cardiovascular failure, Coronary artery disease, Arrhythmias.
- b. Respiratory, chronic obstructive pulmonary disease.
- c. Gastrointestinal, Hepatitis, Malabsorption, Inflammatory bowel disease.
- d. Genitourinary chronic renal failure.
- e. Endocrine, Diabetes, Thyroid disease, Pituitary/adrenal disorders.
- f. Hematological, Anemia, Leukemia, Bleeding clotting disorders.
- g. Musculoskeletal, Arthritis, Osteoporosis.
- h. Neurological, Stroke, Palsy, Mental retardation.



## Risks

Risks associated with the surgical procedure fall into four broad categories:

1. Immediate anesthetic and surgical risks.
2. Psychological and psychiatric risks.
3. Medical threats to long-term retention.
4. Long-term deleterious effects of implants on health.

The risks may include:

Inadvertent perforation of the nasal and maxillary sinus, local and systemic infections, perforation into soft tissue spaces, and nerve injury. Temporary conditions that might result from implant placement may include pain and swelling, speech problems and gingivitis. Long-term problems may include nerve, local or systemic bacterial infections, and



## Important Warning

Lack of adequate practitioner training is a major risk factor in the success of the implant procedure and might endanger patient health. Therefore, no implantation should be performed without prior adequate training from a certified institute.

## Precautions

1. The implant is supplied for one time use. DO NOT RESTERILIZE.
2. Use of implant does not require any unusual preoperative antibiotic prophylaxis.

## Surgical Procedures **Step by Step Protocol**

The surgical manual is designed to provide an overview of the presurgical and surgical procedures applicable to the MIS dental implant system. Successful implantation requires a combination of several factors. This step by step protocol is provided to ensure that the correct surgical procedures are followed.



### **Step 1.**

Patient Selection  
(General medical history)

Patients must be carefully selected for surgery. Medical histories should be evaluated for indications and contra-indications. Medical contra-indications includes: immunodeficiency or immunosuppressive treatment; evolving or active malignancy or following head and neck radiation; poorly controlled diabetes or other hormonal disorders; blood dyscrasias and serious bleeding disorders; recent myocardial infarction, severe cardiac insufficiency and valve pathology; general bone diseases

affecting the jaws; hypersensitivity to specific implant components; psychiatric or personality disorders that limit or interfere with patient understanding and compliance.

### **Caution & Additional Tests**

Caution and additional tests should be used in the following cases: patients with preexisting myocardial disorders, anti-coagulant therapy; arrhythmias and valvular pathology; auto-immune diseases; pregnancy, alcoholism; any respiratory, urinary, gastrointestinal, neurological, hematological, endocrine or other systemic disorders and patients taking bisphosphate medication as a treatment for osteoporosis, mucocutaneous disorder. It is recommended, in cases with questionable medical histories, to consult with the patient's general practitioner and/or appropriate specialist before any surgical procedure.



## Step 2.

### Dental Conditions and Oral Hygiene

A complete and thorough intraoral examination must be performed and recorded. This must include an evaluation of the dentition, oral hygiene, smoking and patient attitude to oral health. Implant procedures should not be performed on patients with active osteolytic, periodontal conditions, active periodontal disease inflammatory or infectious areas at the implant site. Extreme bruxing and clenching forces should be taken into consideration. It should be emphasized that acceptance of a treatment plan that includes dental implants must be preceded by a detailed explanation of what is involved and the influence of heavy smoking (>10/day), uncontrolled diabetes, untreated periodontal disease and bad oral hygiene habits.



## Step 3.

### X-ray Examination (Diagnostic measures)

Panoramic and tomographic radiographs should be obtained. Vertical and horizontal

dimensions of the implant site should be measured and charted. The anatomical relationships of neighboring teeth and proximity to the anatomical structures such as the mandibular canal, maxillary sinus and base of the nose must be evaluated. Bone inclination and shape should also be taken into account. Surgical guide with radioopaque markers is recommended for the tomographic radiographs and for the surgery in complicated cases.



## Step 4.

### Treatment Plan (Patient cooperation)

The various treatment alternatives under consideration should be presented and explained to the patient, i.e. complete removable prostheses retained by ball or bar attachments; screw-retained, fixed or removable crowns or bridges; cemented crowns or bridges, etc. The final choice of treatment must be decided upon with the patient. It is important to obtain patient approval of the treatment plan.

## Surgical Procedures

# Step by Step Protocol



### Step 5A.

#### Surgical Phase

This stage must be done under a strict infection control protocol. Preoperative antibiotics are recommended. Routine blood tests should be performed and blood pressure and pulse checked. emergency resuscitation apparatus should be available. Premedication may be considered for an apprehensive patient or before extensive procedures.



### Step 5B.

#### Implant Selection

Choice of the correct MIS implant should be based on previous evaluations of

bone shape and dimension and take into consideration both functional and aesthetic requirements. The implant is supplied in a sealed and sterilized package. An implant whose sterility has been compromised should not be used. An implant in a partially opened or defective package should not be used. Attention should be paid to the expiration date.

The shape, size and lot number of the implant should be noted (use the attached label). Implant placement should be performed in accordance with the "implant placement procedure" described on pages 20-43.

The sale of this implant is restricted by law to licensed dentists. This surgical procedure should only be performed by someone with appropriate training. Initial planning is of the utmost importance.

At this is a prosthetic driven procedure the dentist performing the prosthetic stage of the treatment should be an active participant together with the surgeon when making decisions affecting the choice of the implant, the type of the prosthesis and the three dimensional positioning of the implant.

# 12 24

## Step 6.

Evaluating Osseointegration

After a post-operative integration period of 12 weeks in the mandible and 24 weeks in the maxilla, the implants may be exposed. Osseointegration is evaluated clinically in conjunction with radiographs.



## Step 7.

Restoration

The final restoration is provided according to the treatment plan, taking into account occlusal and aesthetic requirements. MIS superstructures must be used for MIS implants.



## Step 8.

Follow-up

The follow-up includes an annual check-up including peripical radiograph, re-evaluation and patient education in maintenance and pedantic oral hygiene.



## **Surgical Kits.**

- 56. Surgical Kit Description
- 58. Advanced Surgical Instrument Kit
- 60. Kit Contents

The Surgical Kit  
**Surgical Kit  
 Description**

MIS offers a general surgical instruments kit for two-stage procedures. Ongoing research and development efforts by MIS engineers and scientists have produced this kit to simplify common procedures and provide creative solutions for dealing with complex cases. The surgical kit contains tools for implantation, as well as spare places where drills can be placed according to the doctor's needs.



▲  
 MK-0035, General Surgical Instruments Kit





### Additional advantages of the surgical instrument kit:

- The surgical kit is made of medically approved plastic.
- The surgical kit can be fully sterilized using an autoclave with a temperature that does not exceed 134°C (273°F).
- The surgical kit is small in size, and therefore easy to store.

The modular trays represent the optimal solution in terms of cleaning, decontamination and sterilization due to the absence of hidden surfaces. The steam flow is optimized to maximize the elimination of condensation.



### Warning

Damage to trays:

If the sterilization temperature is greater than 150°C, damage on propylux components is to be expected. Radel, steel and silicone components may support repeated exposures to temperatures up to 180°C, but the lifetime of the trays may be shortened.

When sterilizing cutting edges and the tray, the use of inappropriate chemicals may damage the trays. The trays have to be handled carefully to avoid any breakage. Never use trays with broken areas.



### Cleaning Procedure

It is not recommended to use the following for appliances made of stainless steel:

- Cleaning and disinfection agents containing high rates of chlorine
- Cleaning or disinfection agents containing oxalic acid. It is not recommended to use the following for an appliance or a material marked by a color code:
- Detergents and cleaning agents containing high rates of the aforementioned chemicals
- Extremely high temperature during cleaning and sterilization of the product.

Please Note:

- Please conduct a visual inspection of the appliances. Do not use faulty and dull appliances. Clean and disinfect such appliances separately
- Do not allow traces/ residue (blood, secretion, tissue residue) to dry on the appliances. Always soak in the disinfecting fluid immediately after use
- Use only stainless steel dedicated detergents and strictly follow the usage instructions
- Rinse the disinfectants and cleaning agents thoroughly with water
- Do not store/ put away appliances that are still damp or wet
- Use only nylon bristle brush to remove accumulated dirt residue off the appliance. Clean the cavities and hollow spaces thoroughly
- Clean excessively contaminated appliances using an ultrasonic bathtub
- Do not clean/ disinfect (even in an ultrasonic bathtub) or sterilize together appliances made of different materials
- During mechanical cleaning, please ensure that the parts have no contact with the other parts (to prevent damage)
- After mechanical or manual cleaning, all surgical appliances must be sterilized in an autoclave, at a temperature of 134°C (273°F) at a pressure of ≈315 Kpa, for six (6) minutes. Do not exceed the temperature of 134°C during sterilization. It is not recommended to sterilize in hot air, since you cannot gain adequate control over the temperature
- Do not sterilize rusty appliances
- Inspect for corrosion after sterilization.

The Surgical Kit  
**Advanced  
 Surgical  
 Instrument Kit**

- MK-0035** | Without drills
- MK-IR35** | With internal irrigation drills
- MK-EI35** | With external irrigation drills

MT-PP001

MT-BT020

MT-BT028


MT-BT032

MT-BT038

MT-BT045


MT-BT050

MT-BT055




MT-DE001

MT-DI001



## The Surgical Kit Kit Contents

Surgical Kit include tools that are designed specially for the step by step implantation process. Correct preparation of the implant site ensures efficient and accurate installation. The MIS surgical kit contains the following items:





The Surgical Kit  
**Kit Contents**



			Dimensions	Material
	<b>MT-PP001</b>   Parallel Pin		Ø3.20/2mm length 24mm	Stainless steel
	<b>MT-HW001</b>   Hand Wrench		length 10mm	Stainless steel
	<b>MT-HHK24</b>   Hand Key for Int. Hex. Connection		length 18mm	Stainless steel
	<b>MT-MMH10</b>   Short Motor Adapter For Int. Hex. Connection		length 25.30mm	Stainless steel
	<b>MT-MMH20</b>   Long Motor Adapter For Int. Hex. Connection		length 33.30mm	Stainless steel
	<b>MT-MRH10</b>   Hex. Ratchet Short Adapter Int. Hex. Connection		length 16mm	Stainless steel
	<b>MT-MRH20</b>   Hex. Ratchet Long Adapter Int. Hex. Connection		length 22mm	Stainless steel
	<b>MT-LMK13</b>   Long Motor Key Adapter 0.05"		length 30mm	Stainless steel
	<b>MT-SMK13</b>   Short Motor Key Adapter 0.05"		length 22mm	Stainless steel
	<b>MT-DE001</b>   Drill Extender		length 28.70mm	Stainless steel
	<b>MT-GDN33</b>   Countersink for standard platform implant system		Ø3.75/4.2mm length 26mm	Stainless steel
	<b>MT-GDN50</b>   Countersink for wide platform implant system		Ø5/6mm length 26.8mm	Stainless steel

## The Surgical Kit Kit Contents

			Dimensions	Material
	<b>MT-BT020</b>   Body Try in 2mm		Ø2.00mm length 28.5mm	Stainless steel
	<b>MT-BT028</b>   Body Try in 2.80mm		Ø2.80mm length 28.5mm	Stainless steel
	<b>MT-BT032</b>   Body Try in 3.20mm		Ø3.20mm length 28.5mm	Stainless steel
	<b>MT-BT038</b>   Body Try in 3.80mm		Ø3.80mm length 28.5mm	Stainless steel
	<b>MT-BT045</b>   Body Try in 4.50mm		Ø4.50mm length 28.5mm	Stainless steel
	<b>MT-BT050</b>   Body Try in 5.00mm		Ø5.00mm length 28.5mm	Stainless steel
	<b>MT-BT055</b>   Body Try in 5.50mm		Ø5.50mm length 28.5mm	Stainless steel
	<b>MT-DI001</b>   Implant Direct Indicator		length 15mm	Titanium
	<b>MT-HHR13</b>   Long Hand Screwdriver for 0.05" Hex.		length 22mm	Stainless steel
	<b>MT-HHR07</b>   Short Hand Screwdriver for 0.05" Hex.		length 16mm	Stainless steel
	<b>MT-RI030</b>   Ratchet Wrench		length 75mm	Stainless steel



## Kit Options



**MK-EI35** | Advanced Surgical Kit with external irrigation drills

**MT-TDN19** | Marking drill 1.90mm external irrigation

**MT-TDN20** | Pilot drill 2mm external irrigation

**MT-TDN28** | Twist drill 2.80mm external irrigation

**MT-TDN32** | Twist drill 3.20mm external irrigation

**MT-TDN38** | Twist drill 3.80mm external irrigation

**MT-TDN45** | Twist drill 4.50mm external irrigation

**MT-TDN50** | Twist drill 5mm external irrigation

**MT-TDN55** | Twist drill 5.50mm external irrigation



**MK-IR35** | Advanced Surgical Kit with internal irrigation drills

**MT-TD110** | Marking drill 1.90mm

**MT-TD200** | Pilot drill 2mm

**MT-TD280** | Twist drill 2.80mm

**MT-TD320** | Twist drill 3.20mm

**MT-TD380** | Twist drill 3.80mm

**MT-TD450** | Twist drill 4.50mm

**MT-TD500** | Twist drill 5mm

**MT-TD550** | Twist drill 5.50mm



**MK-0035** | Advanced Surgical kit supplied without drills.





## **Drills.**

- 68. Use of MIS Drills
- 70. Color Code
- 72. Drill Indications
- 75. Preparation of Fixture Site
- 81. Drill Cutting Capability
- 82. Ceramic Drills
- 83. Drill Maintenance

## Drills

# Using MIS Drills

The implant placement procedure includes the use of several drills with different diameters and characteristics. The great variety of drills offered by MIS provide drilling solutions for all the various stages and processes of the implantation site preparation. MIS drills include drills with internal and external irrigation, as well as conical and ceramic drills. Most MIS drills are groove marked for depth control and are color coded for immediate identification of drill diameter. All drilling is carried out under internal irrigation with saline.

### Features

MIS drills are designed to be used with all MIS implants. The drills are available in 3 varieties. In each diameter size there are drills with an internal irrigation hole, without an irrigation hole and a short drill. Each diameter is color coded and a different color is marked on the middle part of the drill. The groove marks on the drills specify the depth of the insertion 6mm, 8mm, 10mm, 11.5mm, 13mm and 16mm. The short twist drills have a 6mm mark included (MT-TDS20-MTTDS55) for our line of 6mm implants. The short drills will have marks at 6mm, 8mm and 10mm. The drills may be used for depths of 11.5mm with the mark being found


at the junction of the color mark and above the 10mm mark. All MIS drills have a 120°C cutting degree. The sharpness and high quality of the drill allows for up to 30 uses. During drilling, the temperature must not exceed 47°C. Each drill has a podium on which the drill stopper can be locked in place.

### Drill Stopper


MIS provides drill stoppers that give the dental surgeon simple and accurate depth control. The sleeves are available for standard implants in 8,10,11.5,13 and 16mm with a different ring color for each diameter (color coded drill matches sleeve color).




Drill stoppers



OUT Proceeding



IN Proceeding




**MT-TD320**  
Twist drill 3.20mm

**MT-TDS20**  
Pilot Short drill 2mm

Drills  
**Color Code**




Color-code is used for easy identification of drills or implants diameters as follows:



## Drills

### Drill Indications

The calculation presented refers to the height between the actual drilling depth and the implant's apex, and shows the difference in length after final drilling. This calculation is important for planning the drilling depth.










- \* The diameter marking drill is used with a continuous in and out motion at the desired location on the alveolar crest, as predetermined with the aid of a radiograph.
- \* Do not drill without irrigation and do not exceed 47°C.
- \* Note: Any procedure recommended by MIS cannot replace the judgment and experience of the surgeon.

## BiOCOM

Calculation of the gap between final drill tip and implant insertion depth:











Biocom Implant	Drills		Gap
 Ø3.30mm	Ø2.80mm		0.2mm
 Ø3.75mm	Ø3.20mm		0.3mm
 Ø4.20mm	Ø3.80mm		0.4mm
 Ø5mm	Ø4.50mm		0.5mm
 Ø6mm	Ø5.50mm		0.6mm

## Drills

### Drill Indications

## SEVEN

Calculation of the gap between the final drill tip and implant insertion depth:

Seven Implant	Drills		Gap
 Ø3.30mm	<sup>1</sup> Ø2.20- <sup>2</sup> Ø3.20mm		0.2mm
 Ø3.75mm	Ø2.80- Ø3.60mm		0.3mm
 Ø4.20mm	Ø3.30- Ø4.10mm		0.5mm
 Ø5mm	Ø4.10- Ø4.90mm		0.6mm
 Ø6mm	Ø5.10- Ø5.90mm		0.95mm

75.

## Drills Preparation of the Fixture Site

### Marking Drill

The Marking Drill is used for creating a reference point in the center of the ridge, to mark the drilling location for further drilling. The Marking Drill is usually inserted to a depth of 1mm without the use of force. Recommended drilling speed is 1200-1500 RPM. The Marking Drill supplied is 34mm in length and 1.90 mm in diameter.

Marking Drill ▶  
MT-TD110



Drills

### Spade Drill

A Spade Drill has a diameter of Ø1.9mm and a sharp tip. Spade Drills are mostly used in Flapless procedures.

Recommended drilling speed: 1200-1500 Rpm. The Spade Drill is 27.5mm in length and made of stainless steel.



◀ Spade Drill MT-SMD10

## Drills

# Preparation of the Fixture Site

### Pilot Drill

The Pilot Drill is the first invasive drill that is used for preparation of the fixture site. The Pilot Drill is used to drill to the desired depth and angle of implantation according to the implant length (8mm,10mm,11.5mm,13mm and 16mm). Recommended drilling speed is 900-1200 RPM.

The Pilot Drill is important, because it determines the initial drilling angle. A parallel pin should be used to check the drilling angle. When placing more than one implant, place a parallel pin in each completed hole before proceeding to the next site. Align the Pilot Drill parallel to the previous pin when available bone permits.

The Pilot Drill is 37mm long and 2 mm in diameter.

**Pilot Drill** ▶  
MT-TD200



### Checking drill depth Implant site depth probe MT-BT110

After the implantation site is drilled, it is recommended to check its depth and suitability to the required implant length. The Implant Site Depth Probe (MT-BT110) is used for this.

\* Made of stainless steel






## Drills

# Preparation of the Fixture Site

### Final Drill

The Final Drill is the last drill used before inserting the implant. The choice of Final Drill is determined by the implant's diameter. The different diameters are: 2.80mm, 3.20mm, 3.50mm, 3.80mm, 4.50mm, 5mm, 5.50mm. The recommended drilling speed is 400-700 Rpm (700 Rpm for small diameter and 400 Rpm for large diameter).

The Final Drills are standard drills used for standard procedures. When drilling the site, they may be used as the last drill. However, in certain cases, there may be an additional instrument used for final preparation of the site. Use of a Final Drill is determined by the specific procedure according to bone and implant type.

	Implant	Length Options	
		Short	Long
 <p>Ø2.80mm</p>	Ø3.30mm	28.4mm	37.5mm
 <p>Ø3.20mm</p>	Ø3.75mm	28.5mm	37.6mm
 <p>Ø3.80mm</p>	Ø4.20mm	28.7mm	38.7mm
 <p>Ø4.50mm</p>	Ø5mm	28.9mm	38.2mm
 <p>Ø5-5.50mm</p>	Ø6mm	28.9mm	38.2mm

## Checking drill depth



After the implantation site is drilled, a check of its depth and suitability to the required implant length is recommended. The Implant Site Depth Probe (MT-BT110) is used for this aim.

Drills

## Preparation of the Fixture Site

### Countersink (MT-GDN33, MT-GDN50)

The Countersink is used to widen the crest to prevent pressure on the implant's neck. Depth marks of 3.75mm, 4.20mm, 5mm and 6mm appear on the Countersink to indicate the suitable implant diameter. The recommended drilling speed is 200-500 RPM.

**Standard**  
MT-GDN33

**Wide**  
MT-GDN50



# SEVEN

## Special Final Drill

A specially designed final drill is recommended for use in bone types 1 and 2 for 6mm, 8mm, 10mm, 11.50mm, 13mm and 16mm SEVEN implants in order to prevent pressure on the implant's neck. The special final drill is supplied with every implant, allowing a short and safe drilling procedure. The recommended drilling speed is 200-400 Rpm.


Final Drill for implant diameters





## Drills

# Drill Cutting Capability



**Test conditions:** Drill diameter: 2 mm  
Drill speed: 600 RPM  
Drill feed: 0.04 mm/rev  
Drill material: fiberglass

Test bench- force transducer: Piezo  
drill feed obtained by DC motor  
controlled by a displacement  
potentiometric transducer

## Conclusion

The MIS stainless steel drills with special sharpening have greater endurance and efficiency in drilling. They minimize the heat generated during the cutting action.

## Drills

# Ceramic Drills

Features of Ceramic drills include reduced vibration, pleasant smooth operation and continuous substance removal, whether used in implantology or oral surgery.

The MIS Ceramic drills are made of a high performance ceramic mix of partly stabilized yttrium and aluminium ceramic. The mixture of these two established materials provides MIS Ceramic drills with an above-average bending strength of 2,000 MPa. In comparison, the bending strength of zirconium oxide ceramic, used among other things in the manufacturing of root posts is 1,200 Mpa.

**Advantages:** Metal-free, biocompatible, corrosion-free



**MT-CRD21**  
Marking Drill

Dimensions:  
Ø2.10mm  
length 28.5mm

Material:  
Zircon-dioxide  
ceramics



**MT-CRD20**  
Pilot Drill

Dimensions:  
Ø2mm  
length 33.5mm

Material:  
Zircon-dioxide  
ceramics



**MT-CRD28**  
Twist Drill

Dimensions:  
Ø2.80mm  
length 35mm

Material:  
Zircon-dioxide  
ceramics

## Drills Drills Maintenance

Correct and careful maintenance of MIS drills is extremely important. Damage to drill tips can cause significant impairment of drill function. The following are detailed instructions for proper maintenance.

### Instructions for Maintenance of Drills Prior to First-Time Surgical Use

**Stage 1:** Light Cleaning and Rinsing - Drills should be dipped in detergent, rinsed, and dried.

**Stage 2:** Sterilization - Drills should be sterilized in an autoclave at 134°C (273°F) at a pressure of 315 Kpa For a 6 minute duration. Do not exceed 134°C during sterilization

**Stage 3:** During Use - Drills should be soaked in a sterile saline solution until the cleaning stage.

### Instructions for Cleaning and Storage of Drills After Use

**Stage 1:** Cleaning - Drills should be brushed with detergent to remove any remaining blood or tissue.

**Stage 2:** Ultrasonic Cleaning - Drills should be cleaned in an ultrasonic bath with appropriate detergent. Note: During ultrasonic cleaning, contact between drills should be avoided.

**Stage 3:** Rinsing - Drills should be rinsed under running water and dried.

**Stage 4:** Lubrication - (Required if more than 4 weeks of storage is expected) Drills should soak for 10 seconds in dental oil, then removed from the solution and left to dry for 30 seconds without rinsing or towel


drying, then placed in the surgical kit.

**Stage 5:** Sterilization - Drills should be sterilized in an autoclave at 134°C (273°F) at a pressure of 315 Kpa For a 6 minute duration. Do not exceed 134°C during sterilization.

**Stage 6:** Storage/Use - At this stage, kits are ready for long-term storage; drills can be used immediately upon opening the kit.

### Recommendations

- Cutting tools should be used for a maximum of 30 drillings.
- Sterilized water should be used in order to avoid surface stains.





# Surgical Tools.

- 86. Specialized Surgical Tools
- 98. Screw Tests
- 99. General Information

Surgical Tools  
**Specialized  
Surgical Tools**

MIS provides a variety of surgical tools to assist the surgeon during the implantation operation. The wide variety of MIS tools are specially designed with the purpose of facilitating a safe, simple and short implantation procedure.



## Torque wrench MT-RI040



### Features

The Torque wrench is designed for tightening or loosening screws and for implant insertion. It also ensures the optimal transmission of force during implant insertion.

The Torque scale ranges from 15-45 Ncm, with an accuracy of plus or minus 5%. The scale on the opposite side can be used as a reverse torque.

The maximal load, as indicated by the scale on the wrench body, should not be exceeded. There is also a scale on the opposite side.

### Advantages


It provides a removal torque check before implant loading.

### User Instructions

1. Connect the torque wrench **A** to the desired key.
2. Place the torque wrench in the mouth.
3. While exerting finger pressure on the handle **B**, turn the torque wrench slowly in a clockwise direction **C** until the desired torque is reached.

### Material

- Stainless steel
- Sterilize before use (in autoclave up to 134°C/273°F)




### Sterilization

- The device is not sterile.
- The device must be sterilized before use by autoclave, at a temperature of 134°C (273°F) and at a pressure of ~315 Kpa for a 6 minute duration. Do not exceed 134°C during sterilization.

Surgical Tools  
**Specialized  
Surgical Tools**

**Implant site depth probe**  
MT-BT110



**Features**

For an easy and accurate measurement of the hole drilled in the bone, this simple procedure helps define the optimal size of implant to use in each case.

Depth of measurement: 8mm, 10mm, 11.5mm, 13mm and 16mm. The depth probe includes

an apical ball that ensures a fast, tactile examination of the cavity when checking implant length.


Dimensions:  $\varnothing 1.80 / \varnothing 2.70\text{mm}$ , length 100mm.



The MIS general surgical kit (MK-0035) contains several unique surgical instruments, as well as tools that can be adjusted for direct connection to the implant.

### Implant Direction Indicator MT-DI001


This surgical instrument reveals the condition of a particular implant by showing the implant direction. The implant indicator is connected directly to the implant and shows the direction of the implant. The implant indicator contains groove marks indicating gingival heights (each groove mark indicates 2 mm of gingival height).



Implant Direction Indicator  
MT-DI001

Surgical Tools  
**Specialized  
 Surgical Tools**

MIS provides a key designed specifically for the extraction of mountless standard or wide implants, placed in very soft bone or in sinus lift procedures. The key can be manipulated manually or with a ratchet. The connection between the key and the implant is facilitated by means of a screw that attaches to the thread of the implant. This allows for a firm for a firm connection between implant and key and for a safe and simple implant extraction.




The key consists of two parts: the body and the screw. The key is composed of stainless steel. It is recommended to dismantle the two key parts (body and screw) before sterilization.

## Direct Hand And Ratchet Hex Key MT-HMR05 / MT-HMR10

Step

1.



### By Hand

Tightening the screw to the implant

Step

2.




### By Ratchet

Ratchet connected to top of the key in order to pull implant

Surgical Tools  
**Specialized  
 Surgical Tools**

**Tools for Int. hex. Connection**

MIS provides two tools that connect directly to the implant hexagon and indicate the hexagonal location of the implant. The upside that connects directly to the internal hexagon is 2.45mm for standard / wide platform and 2.10mm for narrow platform. The insertion tools for implants are also suitable for hand ratchet and motor.



Hand key ▶



MT-HHK24

Ratchet key ▶



MT-MRH10



MT-MRH20



MT-NRH20

Motor key ▶



MT-MMH10



MT-MMH20

Surgical Tools  
**Specialized  
Surgical Tools**

**The Friction fit abutment assembly contents**

MT-IE172/ MT-IE161

The friction fit extractors (MT-IE171 standard/wide and the MT-IE161 narrow) are specifically designed to separate the friction fit abutments from the implant. The extractors are color coded, Blue for standard abutments and Yellow for narrow abutments.



Int. connection  
abutment extractor ▶  
MT-IE172






◀ Int. connection  
abutment extractor, NP  
MT-IE161



## Extractor key

The extractor key is the extractor of friction fit abutments from the implants. Axis force activated on implant axis, take out the abutment from the implant.

### For standard / wide implants



### For narrow implants



Surgical Tools  
**Specialized  
Surgical Tools**

**SOS Broken Screw Kit**  
MT-TF172 / MT-RT001/ MT-HW001

The SOS Broken Screw Kit has been designed to assist the surgeon during the delicate operation of removing a broken screw.



◀ **SOS Broken Screw Kit**  
MK-0041



## SOS Tools.

These tools are recommended for insertion of the implant.

**Thread Form**  
MT-TF172



**Retriever**  
MT-RT001




**Hand Wrench**  
MT-HW001




1.

- A. Connect the retriever to the micromotor.
- B. Adjust the micromotor to low speed (15-25 RPM), max torque and in reverse mode.



2.

- A. Apply pressure with the retriever on the top of the broken screw.
- B. While maintaining the pressure, activate the motor. This action should release the screw. If the screw is still not released, apply intermittent pressure on the screw.



3.


In the case of a damaged thread:

- A. The thread form has to be used carefully.
- B. Be sure to align the thread form parallel to the thread axis.
- C. Thread in the clockwise direction and for every complete turn, release the pressure of the thread form by turning it 30° in a reverse direction and repeat.
- D. In instances where greater torque is in needed, the hand wrench or ratchet may be used.


## Surgical Tools

# Screw Tests

### Tensile test of dental screws



### Fatigue test of dental screw



## Surgical Tools **General Information**

99.

The wide variety of MIS surgical tools requires several maintenance guidelines:



Instrument maintenance for surgical tools:

### **Disinfection**

Immerse instruments immediately after use.  
Use approved agents only.  
Observe manufacturer recommendations regarding concentration/time/material compatibility.

### **Cleaning**

Remove all residues.  
Use an ultrasonic bath.  
Use anticorrosive cleaning agent.  
Thoroughly rinse away cleaning and disinfecting agents with running water.  
Use distilled water to prevent water spots.

### **Drying**

Dry only with: compressed air, hot air, absorbent paper tissue.

### **Examination**

Perform a visual inspection.  
Dispose of damaged instruments.

#### **Check for:**

Breakouts in blades  
Bent instruments  
Corrosion

### **Sterilization**

All dental instruments should be sterilized.  
Use only sterile packages.  
The device is not sterile.  
The device must be sterilized before use by autoclave, at a temperature of 134°C (273°F) at a pressure of ~ 315 Kpa for a 6 minute duration. Do not exceed 134°C during sterilization.

### **Storage**

Store in a dry, dust-proof area.  
Keep instruments separated from chemicals.





**Packaging.**

- 102. Implant Packaging
- 104. Implant Color Code
- 105. Sticker Description
- 106. Implant Mountless Package Handling

## Packaging **Implant Packaging**

MIS puts great effort into providing end users with safe, simple and up-to-date packages. Our products pass very strict tests so that we can provide doctors with the best tools available.



### **The single implant package**

The implant package consists of three sections: an outer carton containing the instructions, an outer plastic tube and inner plastic tube. We recommend that the instructions be read carefully.

Mountless package



## 10 Implants Package

This package contains 10 implants, it can be used as a stand or as a drawer. This packaging is designed for ease of use by the dentist and the distributor.



Insertion of the adaptor



Taking the implant in and out of the tube

## Mountless.


MIS BIOCOM and SEVEN implants are supplied without mounts. These successful implants remain unchanged, but are supplied without the blue mounts, and are packaged in distinctive and easily recognizable boxes. These boxes include a pull tab for quick and easy opening, and a new titanium sleeve, inside the tube, for maximum sterilization. This change improves the handling, speed and safety of the implants. The BIOCOM and SEVEN mountless implants, have their own tools, available as a set (MK-0037) or purchased separately.



Implant

Packaging  
**Implant  
Color Code**

For easy identification of implant diameter, the cap of each outer tube is color-coded as follows:






105.

## Packaging Sticker Description

Between the outer tube and the inner tube there are two stickers. Printed on these stickers are the lot number, the catalog number and the expiration date.



Packaging  
**Implant  
Mountless  
Package  
Handling**

Make sure using physical and visual examination that the implant is of the right type and dimensions for the specific procedure for each patient.





**Fig. 1**

Open the box by pressing on the marked dotted line, and remove the outer tube from the box.



**Fig. 2**

Open the outer tube by pressing down on the lid and turning the tube counter clockwise.



**Fig. 3**

Remove the inner tube from the outer tube. Open the inside tube's cap on the end containing the implant and mount.

Packaging  
**Implant  
Mountless  
Package  
Handling**



**Fig. 4**

Remove the implant label from the outer tube and use for tracking purposes.



**Fig. 5**

Hold the inner tube so that the implant is on the top, and open the cap.



Use one of the following three options to remove the implant from the inner tube:

**Fig. 6A**

A contra-angle hand piece



**Fig. 6B**

A ratchet wrench



**Fig. 6C**

A hand wrench

Packaging  
**Implant  
Package  
Handling**

110.



**Fig. 7**  
Commence implantation  
procedure



**Fig. 8**  
Open the other end of the  
smaller tube. Remove the  
cover screw from the other  
side of the inner tube using  
the MT-HHR13 key



**Fig. 9**

Attach the cover screw to the implant using the MT-HHR13 key



Final Drill



Implant



Cover Screw



Inner Tube with Titanium Sleeve



Outside Tube

## X-Ray Templates

MIS provides X-ray transparencies for measurements and comparisons. This helps the surgeon in the selection of the appropriate implant type, diameter and length according to each individual case.

The following X-ray templates are available:

- |                        |                         |
|------------------------|-------------------------|
| 1- For BIOCOM implants | Cat No. <b>MC-BIOIN</b> |
| 2- For SEVEN implants  | Cat No. <b>MC-SEVIN</b> |














All transparencies include a ruler (100%) in mm and a radiographic implant transparency scale of 100%(real) and 125%.




## Symbols

Key to the symbols on labels and instruction leaflets:

 LOT	Batch code		Manufacturer
 REF	Catalog number		Do not resterilize
	Single use		Do not use if package is damaged
	Attention, see instructions for use		Authorised representative in the European community
	Date of manufacture		By prescription only
	Sterilized using gamma irradiation		Use by








All rights reserved. No part of this publication may be reproduced, transcribed, stored in an electronic retrieval system, translated into any language or computer language, or be transmitted in any form whatsoever, without the prior written consent of the publisher. Warning: Only a licensed dentist should use these products.

MIC-UP207 Rev.3



MIS Implants Technologies Ltd.  
[www.mis-implants.com](http://www.mis-implants.com)

MIS's Quality System complies with international quality standards:  
ISO 13485:2003 - Quality Management System for Medical  
Devices, ISO 9001: 2008 - Quality Management System and  
CE Directive for Medical Devices 93/42/EEC. MIS's products are  
cleared for marketing in the USA and CE approved.